

## Omni Tips

### 7 and 12 mm Disposable Homogenizing Probes



**The Omni Tips™ patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators.**


Proprietary Omni Tips probes are engineered from durable plastics and are so economical they can be disposed of after each use to eliminate cross contamination during homogenization.

The outer tube is made of polycarbonate. The inner shafts are made of Ultem 1000, a polyetherimide. The patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators.

12 mm Omni Tips are suitable for homogenization where sample size or sample hardness dictates that processing with the 7 mm tip is prohibitive.

### Key features




- Economical: <1% of the cost of stainless steel probes
- Eliminates Cross Contamination
- Ideal for PCR, RNA, and DNA sample preparation
- Maintenance free
- Sterile Omni Tips are Electron Beam treated to ensure they are free of contaminants



“With Omni Tips™ I can be certain that all the tissue sample has been homogenized; there is no sample loss. I found that the Omni Tips™ hard tissue version worked as well as the stainless steel probe.”

*Rob T., National Institutes of Health*

**A selection of our best-selling plastic probes**

<b>Omni Tips 7mm Soft Tissue Probes</b>			
		<b>Quantity</b>	<b>Order no.</b>
	Omni Tips plastic homogenizing probes	25 pack	30750
	Omni Tips plastic homogenizing probes	50 pack	32750
	Omni Tips plastic homogenizing probes	100 pack	34750
<b>Omni Tips 7 mm Hard Tissue Probes</b>			
			<b>Order no.</b>
	Omni Tips plastic homogenizing probes	25 pack	30750H
	Omni Tips plastic homogenizing probes	50 pack	32750H
	Omni Tips plastic homogenizing probes	100 pack	34750H
	<b>Sterile</b> Omni Tips plastic homogenizing probes	25 pack	30750H-S
<b>Omni Tips 12 mm Hard Tissue Probes</b>			
			<b>Order no.</b>
	Omni Tips plastic homogenizing probes	25 pack	301250
	Omni Tips plastic homogenizing probes	50 pack	321250
	Omni Tips plastic homogenizing probes	100 pack	341250

**Note:**

A1000SB adapter is required to use Omni Tips on TH motors before 2016 (with black nose ring), Micro and GLH 850.  
 12-120 Omni Tip adapter is required use Omni Tips on TH motors after 2015 (with blue nose ring).

