

Omni Tissue Master

Tissue homogenizer

The Tissue Master 125 is an extremely economical homogenizer that is well suited for customers that are occasionally homogenizing small volume samples with little variation in sample size. The Tissue Master 125 comes complete with a 7 mm stainless steel probe.

The Omni Tissue Master 125 homogenizer complete kit includes 5 mm, 7 mm, and 10 mm stainless steel generator probes, specifically designed for processing volumes of <200 μ L to 100 mL. These probes can handle a wide variety of different sample sizes, making the Tissue Master a great choice for a laboratory with many applications, but a limited budget.


Specifications

- **Processing range:** <200 μ L to 100 mL
- **Speed:** Variable from 500 to 35,000 rpm
- **Power rating:** 125 watt
- **Weight:** 535 g (18.8 oz.)
- **Sound level:** <72 db
- **Warranty:** 1 year motor warranty
- **Standards approval:** CE Approved



Key features

- Powerful 125 watt motor
- Light weight ergonomic design
- Value priced, complete with 7 mm probe
- Complete kit includes 5 mm, 7 mm, and 10 mm stainless steel generator probes




Omni Tissue Master 125 homogenizer complete kit

Includes 5, 7 and 10 mm probes (available individually)

115V - Order number TM125-115

120V - Order number TM125-220



Value priced homogenizing kit

With 5 mm probe - Order number TM125-115-05

With 7 mm probe - Order number TM125-115

With 10 mm probe - Order number TM125-115-10

The Revvity logo, consisting of the word 'revvity' in a lowercase, sans-serif font.